Instrucciones de uso y garantía
Operating instructions and guarantee
Mode d'emploi et garantie
Instruções de utilização e garantia
Manuale d'uso e garanzia
Bedienungsanleitung und Garantie

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902 170 258

Thank you for purchasing our YF63 My Step. All of our products undergo rigorous testing to ensure maximum quality. To use My Step YF63 correctly, first read through the whole manual before using it.

INTRODUCTION

The device My Step YF63 helps you to combat the effects of reduced mobility and physical inactivity, as it produces a smooth movement of the legs that improves blood circulation and activates muscles without leaving home.

GENERAL SAFETY INSTRUCTIONS.

1. Inspect the device before use to ensure that it is in good condition and that the cable is not damaged. If you notice a defect or missing parts, contact the distributor of the product.
2. If you have any health problem, consult your medical doctor before using the device.
3. The device must not be used by children or the disabled without suitable adult supervision.
4. This device is for home use only.
5. If you experience any discomfort while using the unit, cease use immediately and consult your doctor.
6. An excessive use may cause that the unit gets hot. If this happens stop using it and let it get cold.
7. Do not cover the unit during operation (e.g. with a blanket or pillow) because it may overheat and cause damage
8. Do not use the device in damp/humid surroundings and prevent it from getting wet. Do not touch the device with wet hands while it is operating.
9. To avoid injuries and damages keep out your hands and other parts from the moving mechanism when the device is turned on.
10. Do not insert objects inside the unit.
11. Do not stand on the device.
12. Do not grip the appliance by the cable, do not crush it
13. Do not move the device by pulling the cable.
14. To avoid breakage, do not wrap the cord around the device.
15. Switch off the device completely if the appliance is not going to be used for a while and unplug it.
16. Do not unplug the appliance without switching it off first.
17. In the event of a power cut disconnect the plug from the mains supply.
18. To disconnect the device from the main power supply, pull the plug and not the cable.

19. To prevent the risk of fire, damage or electric shock, use the correct voltage
20. Do not expose the device to high temperatures. Keep it in a dry place.
21. If you notice that the device is not operating correctly, do not attempt to repair, dismantle or modify it. Unauthorised repairs may cause injury, in addition to breaking the terms of the warranty. Contact your technical assistance service.

DESCRIPTION

- OFF– OFF position
1 – Speed 1
2 – Speed 2

USE

1. Place the device on a flat surface.
2. Connect the device to the main power supply.
3. To turn on the unit and work on the speed 1 (lower speed) place the key on the position 1.
4. To work on the speed 2 (higher speed) place the key on the position 2.
5. To turn off the unit place the key on the position OFF.

MAINTENANCE AND CLEANING

1. Unplug the device and leave it cool before cleaning it.
2. Clean it with a soft dry cloth. Do not allow the appliance to come into contact with water or any other liquid. Do not immerse the appliance in liquid to clean it. Do not use abrasive cleaners, brushes, gasoline, glass or furniture polishes or solvents to clean it.

TECHNICAL SPECIFICATIONS

MODEL	YF63 MY STEP
VOLTAGE	24V
POTENCIA	20W

Warranty

This product is guaranteed for 2 years in accordance with that stated in Directive 1999/44/EC of the European Parliament and of the Council 25th May, 1999, regarding certain aspects of the sale of consumer goods and associated guarantees, incorporated into Spanish Law by virtue of Law 23/2003, 10th July; Guarantees on the sale of consumer goods.

To this effect, EXERCYCLE, S.L. VAT n° B-01029032, with business address c/ Zurrupitieta 22, Vitoria, Spain, guarantees this product.

Date of purchase.-	Purchased in.-
Purchaser.-	Dealer stamp.-

Please keep the receipt as proof of purchase.

Should you have any queries please do not hesitate to phone our customer service line:

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